

# VICON CAPTURE.U USER GUIDE

## CONTENTS

---

### About Vicon Capture.U 3

### Use the Vicon Capture.U app 4

Requirements for the Capture.U app 5

Download and install the Capture.U app 6

Attach sensors to a subject 7

Choose a capture mode 9

Reset sensors 25

Get help on the Capture.U app 28

Known issues for the Capture.U app 29

### Use Vicon Capture.U Desktop 30

Requirements for Capture.U Desktop 31

Download and install Capture.U Desktop 31

Connect sensors to Capture.U Desktop 32

Download a trial from the sensors 35

© Copyright 2019 Vicon Motion Systems Limited. All rights reserved.

Vicon Motion Systems Limited reserves the right to make changes to information in this document without notice.

Companies, names, and data used in examples are fictitious unless otherwise noted. No part of this publication may be reproduced, stored in a retrieval system, or transmitted in any form or by any means, electronic or mechanical, by photocopying or recording, or otherwise without the prior written permission of Vicon Motion Systems Ltd.

Vicon® is a registered trademark of Oxford Metrics plc.

Vicon Capture.U™ and Vicon Nexus™ are trademarks of Oxford Metrics plc.

IMeasureU is a trademark of IMeasureU Ltd in New Zealand. The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc and any use of such marks by Vicon Motion Systems Ltd is under license. Other product and company names herein may be the trademarks of their respective owners.

IMeasureU Ltd is a wholly owned subsidiary of Vicon Motion Systems Ltd.

For full and up-to-date copyright and trademark acknowledgements, visit

<https://www.vicon.com/vicon/copyright-information>.

Vicon Motion Systems is an Oxford Metrics plc company.

Email: [support@vicon.com](mailto:support@vicon.com) Web: <http://www.vicon.com>



Edit a downloaded trial	48
Manage connected sensors	49
Erase data from sensors	50
Update sensor firmware	51
Get help on Capture.U	53
Known issues for Capture.U Desktop	54

## About Vicon Capture.U

# About Vicon Capture.U


Working seamlessly with Vicon Blue Trident sensors, Vicon Capture.U is a visual application for capturing raw inertial and analytic data, which is not limited by laboratory settings. It offers real-time data overlaid on video, enabling you to analyze movement both on the field and in the lab.

With capture modes for collecting and analyzing raw data, Capture.U enables you to view graphical and statistics information in real time. You can also save video or export data as CSV or PDF files, for further analysis.

For a brief introduction to Capture.U, see the PDF *Vicon Capture.U Quick Start Guide*.

This guide provides information on using both the Capture.U app, for iOS devices; and Capture.U Desktop, for use on Windows PCs and Macs.

- [Use the Vicon Capture.U app \(page 4\)](#)
- [Use Vicon Capture.U Desktop \(page 30\)](#)

 You can also watch [Vicon Capture.U videos](#)<sup>1</sup>, covering all aspects of using Capture.U, on YouTube.

---

<sup>1</sup> <https://youtu.be/kRuio13OqQk>

Use the Vicon Capture.U app

## Use the Vicon Capture.U app

For information on installing and using the Capture.U app, see these topics:

- [Requirements for the Capture.U app \(page 5\)](#)
- [Download and install the Capture.U app \(page 6\)](#)
- [Attach sensors to a subject \(page 7\)](#)
- [Choose a capture mode \(page 9\)](#)
- [Reset sensors \(page 25\)](#)
- [Get help on the Capture.U app \(page 28\)](#)
- [Known issues for the Capture.U app \(page 29\)](#)

## Use the Vicon Capture.U app

### Requirements for the Capture.U app

The Vicon Capture.U app is available for iOS devices only.

Before you begin, ensure your iOS device meets or exceeds these requirements:

- iOS 10 and BLE 4.2 devices:
  - iPad Air 2, iPad 2017, iPad 2018, iPad Mini 4 and early iPad Pro
  - iPhone 6 and iPhone 7

For optimum performance, this specification is recommended:

- iOS 12 and BLE 5.0 devices:
  - iPad Air 2019, iPad mini 2019, iPad Pro 11-inch, iPad Pro 12.9
  - iPhone 8 and iPhone X and later

This table gives full details of compatible iOS devices:


iOS devices (iOS 10-12)	Bluetooth 4.2	Bluetooth 5
iphones	iPhone 6 iPhone 6 Plus iPhone 6s iPhone 6s Plus iPhone 7 iPhone 7 Plus	iPhone Xs Max iPhone Xs iPhone Xr iPhone X iPhone 8 iPhone 8 Plus
iPads	12.9-inch iPad Pro 2nd generation 12.9-inch iPad Pro 1st generation 10.5-inch iPad Pro 9.7-inch iPad Pro 9.7-inch iPad 6th generation (2018) iPad 5th generation (2017) iPad Air 2 iPad mini 4	12.9-inch iPad Pro 3rd generation 11-inch iPad Pro iPad Air 3rd generation iPad mini 5th generation

Use the Vicon Capture.U app

## Download and install the Capture.U app

To install Capture.U:

- Download the Vicon Capture.U app  from the App Store to your iOS device.

 Capture.U supports only Vicon Blue Trident sensors. Blue Thunder sensors are not supported.

Use the Vicon Capture.U app

## Attach sensors to a subject

The method you use to attach the sensors to your subject depends on the type of movement you want to capture.

- To capture the movement of a subject's limbs, you can attach the sensors using the supplied straps.
- To capture the movement of other parts of the anatomy (for example, where movement of the vertebrae is of interest), you can attach sensors directly, using suitable tape.

**To attach sensors to a subject using the supplied straps:**

1. Insert each sensor into its strap with the IMU symbol facing outwards.



2. Attach the straps to the subject, ensuring that:
  - a. The strap sits snugly against the limb.
  - b. The sensor is oriented with the top (head) pointing in the same direction as the movement of the subject.
  - c. The flashing LED is at the top of the strap.


The following example shows a strap attached so that the sensor sits directly on the medial aspect of the tibia, just above the medial malleolus:



## Use the Vicon Capture.U app

To attach sensors to a subject where straps are not used:

1. Position the sensor on the subject, depending on the required sensor orientation.




2. Secure the sensors using your preferred type of tape (for example, hypoallergenic, double-sided, micropore surgical tape).



Use the Vicon Capture.U app

## Choose a capture mode

On your iOS device, start Capture.U and tap to choose a capture mode.



Your choice of capture mode depends on the number of sensors you want to use, the range you require and the type of output you want:


- For unlimited range and up to 14 sensors (depending on device), with reference video data directly from the device, choose To Sensor. After you have captured the data, download it to your computer with Capture.U Desktop for further analysis.
- For data that is saved directly to your device (where the capture range is constrained by the Bluetooth® range of the device), choose one of the other capture modes.

The capture modes are explained in more detail in the following sections:

- [Capture to sensor \(page 10\)](#)
- [Capture to device \(page 16\)](#)
- [Capture in Real-Time Insight mode \(page 18\)](#)
- [Run Activity Widget demo \(page 24\)](#)

## Use the Vicon Capture.U app

### Capture to sensor



Choose the **To Sensor** mode when you want a high sensor count and unlimited range. After capturing the required movement, you can download your data with Vicon Capture.U Desktop.

An example of this type of usage might be to monitor the track performance of multiple athletes who are running a marathon. After capturing data from all the athletes, you can download the data from the sensors to your computer with Capture.U Desktop (see [Use Vicon Capture.U Desktop \(page 30\)](#)).

 Watch the Vicon video, [Capture to Sensor<sup>2</sup>](#) on YouTube.

Summary of **To sensor** mode:


<b>No. of sensors:</b>	Up to 14, depending on device
<b>Data type:</b>	IMU data
<b>Video:</b>	Reference video
<b>Axes:</b>	12 High G and Low G accelerometer Gyroscope and Magnetometer (option to capture all axes)
<b>Output:</b>	Data (.csv, with Desktop) Video (.mp4) and notes (.txt)
<b>Range:</b>	Unlimited
<b>Capture rate:</b>	High G accelerometer (default): 1600 Hz Low G accelerometer and Gyroscope: 1125 Hz Magnetometer: 100 Hz

<sup>2</sup> [https://youtu.be/yR6tCPH\\_iJU](https://youtu.be/yR6tCPH_iJU)


## Use the Vicon Capture.U app

### To capture to sensors:

1. Open the Capture.U app.
2. In the Capture section tap **To Sensor**.
3. Either create a new session or select an existing one:
  - To create a new session:
    - i. At the top of the screen, tap the plus sign + to the right of **Sessions**.




- ii. In the **Session Name** field, enter a name for your new session.
- iii. If you want to add an image as an identifier for the session, tap **Take Photo** and add the required image.
- iv. Tap **Done**.  
Your new session is displayed in the list of sessions on the left.  
(To make any changes, slide the session towards the left and tap **Edit**.)
- To select an existing session:
  - Tap an existing session in the list on the left.




## Use the Vicon Capture.U app

- At the bottom right, tap **New Trial**, then enter a name for your trial and from the options below, select the data to stream (**Accelerometer ( High G or Low G )**, **Gyroscope**, **Magnetometer**), and whether you want to capture a video.



If you turn on **Capture Video**, you can choose to record **High** (1920 x 1080), **Medium** (480 x 360, the default), or **Low** quality (192 x 144) video.

- When you have made your selections, check below to make sure the data mode is correct (in this case, it is **Data captured to sensor**), and tap **Start Capture** at the bottom right of the screen.



## Use the Vicon Capture.U app

6. Select the sensors, by doing one of the following:
  - Tap to select the required sensors from the list on the right.  
or
  - Tap **Enable Tap to Select** and then tap the required sensor(s) twice.




**i** If one or more of the required sensors is displayed in orange and is labeled **Reset required**, see [Reset sensors \(page 25\)](#).

7. At the bottom right, tap **Start Capture**.  
If you chose to record a video, it is displayed at the bottom left of the screen.  
If Video is not selected, you can capture two trials simultaneously.  
To add notes during capture, tap **Notes** (top right).

## Use the Vicon Capture.U app

8. When you have captured the required data, tap **Stop Capture**.  
Your new trial is displayed in the list of trials on the right of the screen (the most recent trial is at the top).




If you recorded a video, you can play it back on your iOS device without exiting Capture.U by tapping the video icon on the right.



## Use the Vicon Capture.U app

To add notes after you have finished capturing, slide the trial towards the left and then tap **Notes**.



9. To export any video and/or notes that you created, slide the trial towards the left, and then tap **Export**.  
If you recorded a video, you are asked if you want to include it in the export. Tap **Yes**, then select the required export option.  
A .zip file with the name of your trial, which contains any notes (.txt) and video (.mp4) is exported.

✓ Use Quicktime or another suitable video player, eg VLC Media Player, to view the downloaded video file.

10. Download your trial data (.csv) with Capture.U Desktop (see [Use Vicon Capture.U Desktop](#) (page 30)).

## Use the Vicon Capture.U app

### Capture to device



Choose the **To Device** mode when you need a higher sensor count than **Real-Time Insight** mode (which is restricted to two sensors). Data is captured within the Bluetooth range of the device only. In this mode, you can also use reference video.

 Watch the Vicon video, [Capture to Device](#)<sup>3</sup> on YouTube.

Summary of **To device** mode:

<b>No. of sensors:</b>	Up to 6
<b>Data type:</b>	IMU data
<b>Video:</b>	Reference video
<b>Axes:</b>	12 High G or Low G Accelerometer, Gyroscope and Magnetometer
<b>Output:</b>	Data (.csv), video (.mp4) and notes (.txt)
<b>Range:</b>	Bluetooth: 20 m+
<b>Capture rate:</b>	High G accelerometer (default): 800 Hz Low G accelerometer and Gyroscope: 562 Hz Magnetometer: 100 Hz

#### To capture to a device:


1. Open the Capture.U app.
2. In the Capture section tap **To Device**.
3. Add a new **Session** or select an existing one (see [Capture to sensor \(page 10\)](#), Step 3).
4. Create a **New Trial** (see [Capture to sensor \(page 10\)](#), Step 4, but note that the capture options differ, as shown in the above table) and then tap **Start Capture** at the bottom right of the screen
5. Select the required sensors (see [Capture to sensor \(page 10\)](#), Step 6), and then tap **Start Capture**. To add notes during capture, tap **Notes** (top right).
6. When you have captured the required data, tap **Stop Capture**.  
Your new trial is displayed in the list of trials on the right of the screen (the most recent trial is at the top).

<sup>3</sup> <https://youtu.be/uFW21vz9CZE>



## Use the Vicon Capture.U app


Notice that in the filter at the top, **Device** is automatically selected, to display trials recorded directly to your iOS device. If required, you can change this to display trials captured to sensors, or all trials.



You can play back a recorded video on your iOS device without exiting Capture.U by tapping the video icon in the list of trials.



To add notes after capturing, slide the trial left and tap **Notes**.



7. To export your data and any video/notes, slide the trial left and tap **Export**. If you recorded a video, you are asked if you want to include it in the export. Tap **Yes**, then tap the required export option. A .zip file containing your data (.csv), notes (.txt) and video (.mp4) is exported.

✓ Use Quicktime or another suitable video player, eg VLC Media Player, to view the downloaded video file.

Use the Vicon Capture.U app

## Capture in Real-Time Insight mode



Choose **Real-Time Insight** mode to capture movement for real-time data streaming from two sensors with video overlay. You can also define your own benchmarks (that is, you can set a threshold for alerts). For example, you can choose to be notified by a sound when a specified level is reached, so that you can easily recognize when an athlete hits a ball hard enough.

An example of using Real-Time Insight mode might be a coach and an athlete who want to review performance data together over multiple trials on an iOS device, then have the coach send the athlete a report of what was reviewed.

 Watch the Vicon video, [Using Real-Time Insight<sup>4</sup>](https://youtu.be/zPk-WIS5OmY) on YouTube.

Summary of **Real-Time Insight** mode:


<b>No. of sensors:</b>	2
<b>Data type:</b>	IMU data
<b>Video:</b>	Video data real-time overlay
<b>Axes:</b>	3 High G Accelerometer (default) or Low G acceleration or Gyroscope
<b>Output:</b>	Data (.csv), video (.mp4) and report (.pdf)
<b>Range:</b>	Bluetooth: 20 m+
<b>Capture rate:</b>	High G accelerometer (default): 800 Hz Low G accelerometer and Gyroscope: 562 Hz

<sup>4</sup> <https://youtu.be/zPk-WIS5OmY>

## Use the Vicon Capture.U app

To use Real-Time Insight mode:


1. Open the Capture.U app.
2. In the Capture section tap **Real-Time Insight**.
3. At the bottom left of the screen, tap **Select Sensors**.



## Use the Vicon Capture.U app


4. Select up to two sensors (for details, see [Capture to sensor \(page 10\)](#), Step 6), and then tap **Set**.





5. Select either **Acceleration** (x, y, z, or **Resultant** and **Low G** or **High G**) or **Rotation** (x, y, or z).  
Note that **Resultant** is the peak resultant acceleration, which is calculated by  $\sqrt{x^2 + y^2 + z^2}$   
With the **High G** accelerometer selected, the resultant acceleration can reach up to 346 G.  
With the the **Low G** accelerometer selected, the resultant acceleration can reach up to 26 G.
6. At the bottom right of the screen, tap **Start**.
7. To play a sound when a specified level is reached, in the **Benchmarks** area, enter the required values in the benchmark fields.  
The values are indicated by a shaded area on the graph.
8. To view and optionally record real-time video (recorded at 1920 x 1080), tap the Video icon  at the top right of the screen.  
The graph is overlaid on the video. To make the graph easier to see, the

## Use the Vicon Capture.U app

controls at the top right of the screen, enable you to change the opacity of the graph or select a dark or light version.



9. To see statistics for the capture in real time, tap the Statistics icon , towards the bottom right of the screen.  
The minimum and maximum values detected by the sensors are displayed, in addition to how many times the benchmark values have been reached.





10. When you have captured the required movement, tap **Stop**.

## Use the Vicon Capture.U app


11. From the export icons, select the required export type (file or video):

- To save a video to your device, tap the Video icon (top right) and then tap **Save**.



Find saved videos on the **Photos** tab in the Photos app .


- To export either a report (.pdf) or data (.csv), or to browse previously recorded videos, tap the File icon (top right).



## Use the Vicon Capture.U app

- If you choose **Export to PDF**, you can add a title in the PDF preview that is displayed.

The preview shows the information that will be displayed in the exported PDF, with statistics information at the top of the page, a Focused view of the graph (showing the region of interest) and Full Range view (showing the whole range). This lets you check that the graph is as required before you export it. If it isn't, tap **Cancel**, change the **Acceleration** and/or **Rotation** to show a graph of the required view (X, Y, Z, or Resultant) and tap **Export to PDF** again. When you're happy with the preview, tap **Export** at the top to export the PDF.



- If you choose **Export to CSV**, the data values exported can be X, Y and Z, but not Resultant.
- If you choose **Browse Video**, previously captured Real-Time Insight videos are listed. Tap to select a video, then play a minimized preview by tapping the Play icon at the bottom, or tap **Use** at the top right to play in full-screen mode.

## Use the Vicon Capture.U app

### Run Activity Widget demo

The Activity Widgets provide pre-defined workflows that have output metrics.

Activity Widgets might be particularly useful in situations where a quick assessment with minimal/no user setup is required, for example, in undergraduate research projects, grass-roots sports, etc.

For a quick demonstration, you can view an animated GIF example. To do this:

1. From the **Activity Widgets**, select **Swimming**.
2. Hold the iOS device in landscape orientation.
3. Tap the screen to run the example through each workflow screen.
4. When you have tapped through all the screens, you are returned to the home screen.

This provides an overview of the workflow for the activity widgets.



To exit the Swimming GIF example at any time, tap **Home**.  
You are returned to the Capture.U app home screen.


Further, fully interactive Activity Widgets are planned.



Use the Vicon Capture.U app

## Reset sensors

When you select a capture mode and tap **Select Sensors**, if one or more of the sensors needs a reset, it is displayed in orange and is labeled **Reset required**.



This may happen, for example, if you close the app before stopping a capture.

You can reset sensors using the **Reset sensors** option.


To reset sensors:

1. Tap the Settings icon (the cog symbol, bottom right).



## Use the Vicon Capture.U app

### 2. Tap Reset sensors.




### 3. Select the sensor by doing one of the following:

- Select **Tap to Select** and then tap the sensor twice.  
or
- In the list in the app, tap the sensor name.

The sensor is selected.

## Use the Vicon Capture.U app

4. Tap **Reset Selected**.



The reset sensor(s) are now available and the app returns to the **Settings** window.

5. To close the **Settings** window, tap the **Back** button or tap a blank area of the Capture.U screen.

Use the Vicon Capture.U app

## Get help on the Capture.U app

1. Tap the Settings icon (the cog symbol, bottom right).



2. Tap **Help**.
3. Tap one of the options:
  - **Customer Support** Opens the [Vicon Support webpage](https://www.vicon.com/support)<sup>5</sup>.
  - **Privacy Policy** Displays the [Vicon Capture.U legal policy](http://legal.vicon.com/app_privacy_policy.html)<sup>6</sup>.
  - **Terms of Use** Displays the [Vicon Capture.U terms of use](http://legal.vicon.com/ios_terms_of_use.html)<sup>7</sup>.
4. To close the **Settings** window, tap the **Back** button or tap a blank area of the Capture.U screen.

<sup>5</sup> <https://www.vicon.com/support>

<sup>6</sup> [http://legal.vicon.com/app\\_privacy\\_policy.html](http://legal.vicon.com/app_privacy_policy.html)

<sup>7</sup> [http://legal.vicon.com/ios\\_terms\\_of\\_use.html](http://legal.vicon.com/ios_terms_of_use.html)


Use the Vicon Capture.U app

## Known issues for the Capture.U app

Issue	Workaround
When Capture.U is used on an iPhone, it remains in portrait orientation, even when the phone is rotated.	None, as this is intended behavior. When Capture.U is used on iPhones, for optimal use of the screen proportions, it is displayed in portrait mode only.
On iPhones, graphs are always exported in portrait orientation (side-by-side).	None, as this is intended behavior. As phones are locked to portrait orientation (see above), graphs are always exported in portrait orientation (side-by-side).
In Real-Time Insight mode, when you set benchmarks, you can't close the Benchmark keyboard.	To close the Benchmark keyboard, tap the Real-Time Insight graph.
If trial and/or session names are followed by a double space, a period is automatically added. The period is an illegal character, so if you try to use a trial or session name that is followed by a period, the Capture.U app does not save it.	Do not append a double space to trial or session names.
Downloaded High G accelerometer data may exhibit "latching". This is where, when the data becomes close to zero, its value is automatically rounded to zero.	None

## Use Vicon Capture.U Desktop

# Use Vicon Capture.U Desktop



Vicon Capture.U Desktop lets you download data that you collected with Vicon Blue Trident IMU sensors and access a walk-through and documentation. You can export data as video, CSV or PDF files for further analysis.

 Watch the Vicon video, [Using Capture.U Desktop<sup>8</sup>](#) on YouTube.

When using Capture.U Desktop, note that:

- Capture.U Desktop is a fixed size window so can't be minimized or maximized.
- To close Capture.U Desktop, click the **Exit** button (bottom left of the window).

For more information about Capture.U Desktop, see:

- [Requirements for Capture.U Desktop \(page 31\)](#)
- [Download and install Capture.U Desktop \(page 31\)](#)
- [Connect sensors to Capture.U Desktop \(page 32\)](#)
- [Download a trial from the sensors \(page 35\)](#)
- [Edit a downloaded trial \(page 48\)](#)
- [Manage connected sensors \(page 49\)](#)
- [Erase data from sensors \(page 50\)](#)
- [Update sensor firmware \(page 51\)](#)
- [Get help on Capture.U \(page 53\)](#)
- [Known issues for Capture.U Desktop \(page 54\)](#)

<sup>8</sup> [https://youtu.be/\\_RC28CXasrk](https://youtu.be/_RC28CXasrk)

## Use Vicon Capture.U Desktop

### Requirements for Capture.U Desktop

The following requirements are the minimum that are recommended and fully supported:

#### Capture.U Desktop for Windows

- Windows 10 (64 bit)
- NET Framework 4.5

#### Capture.U Desktop for OSX

- MacOS Sierra (10.12)

### Download and install Capture.U Desktop

To install Capture.U Desktop on your computer:

1. Check that your computer software meets the recommended specification (see [Requirements for Capture.U Desktop \(page 31\)](#)).
2. Download Capture.U Desktop from [www.vicon.com/downloads/software](http://www.vicon.com/downloads/software)<sup>9</sup>.
3. Unzip and double-click the .exe file to install.

---

<sup>9</sup> <http://www.vicon.com/downloads/software>

## Use Vicon Capture.U Desktop

### Connect sensors to Capture.U Desktop


To connect sensors:

1. Insert each Blue Trident sensor into its USB adapter.
2. With the supplied micro-USB cables, connect the adapter(s) to the PC or Mac.



If you can't connect all the sensors that were used for capturing data because you used more sensors than you have USB ports, connect the sensors to the available ports. You will be able to connect the remaining sensors later. For more information, see [Download a trial when some sensors are not connected](#) (page 42).


3. On your computer, start Capture.U Desktop.





## Use Vicon Capture.U Desktop





4. To confirm that one or more sensor(s) is connected to Capture.U Desktop, click either **Process** or **Sensors**.  
All connected sensors are displayed.



## Use Vicon Capture.U Desktop

### Capture.U icons on the Captures tab

The following icons are displayed on the **Captures** tab to indicate the status of the sensors:

Icon	Meaning
 IMU-186	Blue Trident sensor is connected to Capture.U Desktop.
 IMU-4	Blue Trident sensor is not connected to Capture.U Desktop.
 IMU-3	Blue Trident sensor was selected for the trial, but no data has been recorded or data has been erased.
 IMU-31	Blue Trident sensor data has been downloaded for the trial.

## Use Vicon Capture.U Desktop


### Download a trial from the sensors

The following steps describe how to download a trial from the sensors where all the sensors used for the capture are connected and no data is missing. For information on other scenarios, see also:

- [Download a trial when some sensors are not connected \(page 42\)](#)
- [Download a trial when sensor data has been erased \(page 46\)](#)

To download a trial from the sensors:

1. Ensure your Blue Trident sensors are connected to your computer (see [Connect sensors to Capture.U Desktop \(page 32\)](#)).  
If a red circle is displayed in the top right of the **Sensors** symbol, the Blue Trident sensors need a firmware update.




Before attempting to continue, update the firmware (see [Update sensor firmware \(page 51\)](#)).


Note that you can't download trials that were captured with the **To Sensor** capture mode until the sensor firmware has been updated.

## Use Vicon Capture.U Desktop

2. Select **Process** and on the **Captures** tab, under **Select trial to download**, click the required trial.




The **Process** screen changes from the **Captures** tab to the **Preview** tab.




## Use Vicon Capture.U Desktop

3. If required, before downloading the trial, you can edit the capture to rename it, crop it, add notes, and/or change the Save location:
  - To rename the trial, click in the **Session** and/or **Trial** fields.




- To crop the trial, so that you download only the region of interest:
  - i. Select the **Crop** button at the bottom right of the graph.



- ii. Select the required part of the trial by clicking on the graph then dragging over the region of interest.  
The selected region is indicated by a rectangle. If you need to adjust the size of the selected region, drag the sides of the rectangle.

## Use Vicon Capture.U Desktop

In the following example, the selected region excludes the start and end of the trial.




iii. Select **Apply** to crop the trial.




## Use Vicon Capture.U Desktop

The time of the trial changes to the cropped time.




- To add notes, click in the **Notes** field (bottom left).




## Use Vicon Capture.U Desktop

- To change the location for downloaded trials, click **Save to**. This displays the File Explorer, enabling you to select the required folder.



4. When you have finished editing the trial, click **Download** to save the trial data to your computer.




An orange progress bar indicates the download status.




## Use Vicon Capture.U Desktop

When the trial has finished downloading, the trial is labeled **Downloaded** and green check marks are displayed on each sensor.



The trial data and any notes you created are saved in the **Save to** folder.




## Use Vicon Capture.U Desktop

### Download a trial when some sensors are not connected

If you can't connect all the sensors that were used for capturing data because you used more sensors than you have USB ports, connect the sensors to the available USB ports. You will be able to connect the remaining sensors later.

To download a trial when some sensors are not connected:

1. Connect sensors to the available USB ports on your computer (see [Connect sensors to Capture.U Desktop \(page 32\)](#)).
2. Select **Process** and on the **Captures** tab, notice that connected sensors are displayed in blue and sensors that are not connected are displayed in gray. This example shows a trial that was captured with four sensors, but with only two of the sensors connected.




3. Under **Select trial to download**, click the required trial.

## Use Vicon Capture.U Desktop

On the **Preview** tab, a message warns you that not all sensors are connected.




4. Edit the trial if required (see Step 3 in [Download a trial from the sensors \(page 35\)](#)). To ensure data formatting is consistent, any cropping of the data of the connected sensors is automatically applied when you connect the remaining sensors.
5. When you have finished any editing, click **Download**.  
When the trial has been downloaded, it is labeled **Partially Downloaded**.



## Use Vicon Capture.U Desktop

6. Connect the remaining sensors to the computer.



7. Select the trial.

On the **Preview** tab, notice that the **Crop** button is unavailable, to ensure the data is formatted in the same way as the previously connected sensors.


If you cropped the trial when you connected the first set of sensors, the trial is automatically cropped in the same way as the previously connected sensors.



8. Click **Download**.

## Use Vicon Capture.U Desktop


When data from all the sensors has been downloaded for the associated trial, green check marks are displayed for each sensor.



## Use Vicon Capture.U Desktop

### Download a trial when sensor data has been erased

1. If you deleted data from one or more sensors that were used in a trial before downloading the data, when you connect the relevant sensor(s), any sensor with deleted data is displayed with a red cross on the **Captures** tab:




You can still download data from the other sensors, but when you select the trial, a message will warn you that some data is missing and ask if you want to download anyway.

2. To continue, select **Download anyway**.


## Use Vicon Capture.U Desktop

3. Edit the trial if required (see Step 3 in [Download a trial from the sensors \(page 35\)](#)) and then click **Download**.

When the trial has been downloaded for all sensors, on the **Preview** tab, you are warned that the trial has not been discovered on all the required sensors.



On the **Captures** tab, the trial is labeled **Partially Downloaded**.




## Use Vicon Capture.U Desktop

### Edit a downloaded trial

To change a trial that you have downloaded (for example, if you want to crop the trial), provided that you haven't deleted the data from the sensor, you can re-download the trial and make the required changes.

1. With the relevant sensor connected, select **Process** and on the **Captures** tab, select the required trial.

A message notifies you that this session has already been downloaded and asks if you want to continue.



2. Click **Continue**.

You can then edit the trial as required and download in the same way as the original download (see [Download a trial from the sensors \(page 35\)](#)).

**i** If you haven't changed the session/trial name and Save location, the trial export files are only differentiated by their download time.




## Use Vicon Capture.U Desktop

### Manage connected sensors

To manage connected sensors, in the main Capture.U Desktop window, click **Sensors**.

From here you can [update sensor firmware](#) (page 51), [erase data from sensors](#) (page 50), and display detailed information about each one.

At the top right of the window, the number of connected sensors is displayed.



For each connected sensor, the following information is displayed (left to right):


- Sensor ID number
- Percentage of storage available for the sensor
- Battery percentage. Connected sensors are charged if they are not already fully charged.

For further information, click the downward arrow.



The following details are displayed:

- **Sessions:** Number of recorded trials on the sensor
- **Firmware version**




## Use Vicon Capture.U Desktop

### Erase data from sensors


To remove unwanted data from your sensors, you delete it using the **Erase** button in the **Sensors** screen.

To erase data from a sensor:

1. Ensure the sensor(s) from which you want to erase data is connected.
2. Click **Sensors** and then for the required sensor, click **Erase**.  
A message notifies you that all data will be erased from the sensor and you're given the opportunity to confirm or cancel.



3. To delete the data from the sensor, click **Erase**.  
When the data has been erased from the sensor, the **Erase** button is disabled and if you click the down arrow, you can see that the **Sessions** information now shows a zero.




On the **Captures** tab, a red cross is displayed on the connected sensor



## Use Vicon Capture.U Desktop


### Update sensor firmware

When Blue Trident sensors need a firmware update, in the main Capture.U Desktop window, a red circle is displayed in the top right of the **Sensors** symbol.




To update sensor firmware:

1. Click **Sensors** and then click **Upgrade firmware**.




## Use Vicon Capture.U Desktop

The status changes to **Firmware upgrade in progress**.



When the firmware upgrade is complete, the sensor details are displayed.




2. If other sensors need to have their firmware upgraded, repeat the above step for each one.
3. You can now download trials that were captured using the **To Sensor** capture mode (see [Download a trial from the sensors](#) (page 35) and [Capture to sensor](#) (page 10)).

## Use Vicon Capture.U Desktop

### Get help on Capture.U

To get help on using Capture.U:

1. Start Capture.U Desktop and click **Help**.



2. Choose from the following options:

- **Guides\*** Visit the [Vicon IMU documentation web page](https://docs.vicon.com/display/IMU/IMU+documentation)<sup>10</sup>.
- **FAQs\*** Visit the [Vicon FAQs web page](https://www.vicon.com/faqs/)<sup>11</sup>.
- **Support\*** Visit the [Vicon Support web page](https://www.vicon.com/support/)<sup>12</sup>.
- **Start walkthrough.** Launch a quick tour of Capture.U Desktop.  
To use this option, ensure that you have first captured data using the **To Sensor** capture mode in the Capture.U app and connected the sensor(s) to the PC.  
To end the walkthrough at any stage, click **Help** and then click **End walkthrough**.
- **Version check.** Click the version number (bottom right) to display more detailed version information about Capture.U Desktop.

\* Internet access required

<sup>10</sup> <https://docs.vicon.com/display/IMU/IMU+documentation>

<sup>11</sup> <https://www.vicon.com/faqs/>

<sup>12</sup> <https://www.vicon.com/support/>

## Use Vicon Capture.U Desktop

### Known issues for Capture.U Desktop

Issue	Workaround
When installed on a clean Windows machine, (ie, on a machine on which you haven't previously installed any Vicon software), Capture.U Desktop starts up with an error and shuts down.	Install <code>vc_redist.x64.exe</code> , which is included in the Capture.U Windows installer.
Downloaded High G accelerometer data may exhibit "latching". This is where, when the data becomes close to zero, its value is automatically rounded to zero.	None