
What's New in Vicon Polygon 4.3.1

Contents

About this guide	2
Vicon Polygon product description	3
Requirements for Vicon Polygon	4
License server requirements	5
New features of Vicon Polygon	6
Video title display	7
Time Bar displays active data folder	8
Other updates	9
Vicon Polygon addressed issues	10
Vicon Polygon known issues	11
Contact Vicon	12

© Copyright 2016 Vicon Motion Systems Limited. All rights reserved.

Vicon Motion Systems Limited reserves the right to make changes to information in this document without notice.

Companies, names, and data used in examples are fictitious unless otherwise noted. No part of this publication may be reproduced, stored in a retrieval system, or transmitted in any form or by any means, electronic or mechanical, by photocopying or recording, or otherwise without the prior written permission of Vicon Motion Systems Ltd.

Vicon® is a registered trademark of OMG plc. Vicon Blade™, Vicon Control™, Vicon Lock™, Vicon Lock+™, Vicon Nexus™, Vicon MX™, Vicon Polygon™, Vicon Studio™, T-Series™, Vicon Tracker™, Vicon Vantage™, Vicon Vero™ and Vicon Vue™ are trademarks of OMG plc.

VESA® is a registered trademark owned by VESA (www.vesa.org/about-vesa/). Other product and company names herein may be the trademarks of their respective owners.

Vicon Motion Systems is an OMG plc company. Email: support@vicon.com Web: <http://www.vicon.com>

About this guide

This document includes the following information about Polygon 4.3.1.

- Product information that you need to get started using this release of Polygon
- New features and functionality included in this release
- Issues fixed in this release

Vicon Polygon product description

Polygon offers unique intuitive and time-saving tools for reporting and presentation, to enable you to communicate and share data.

Polygon 4.3.1 is a point release that includes new features for enhanced usability. For details see [New features of Vicon Polygon on page 6](#).

Requirements for Vicon Polygon

Vicon recommends the Windows 7 (64-bit) operating system for use with Polygon 4.3.1. Polygon has also undergone limited testing under the Windows 10 (64-bit) operating system. Although Polygon may install and function under other Microsoft Windows operating systems, this is not officially supported or recommended by Vicon.

Polygon is fully integrated with Vicon Nexus and Vicon BodyBuilder through Eclipse data management.

License server requirements

If your organization licenses its Vicon product(s) from a server, ensure that the server meets the following criteria:

- It is unlikely to be powered down.
- It has a static IP address, so that client machines can find it easily and reliably.
- It is unlikely to have its hardware configuration changed, especially the HDD and NICs.

 **Tip**


If your organization already uses a server, it is likely to fulfill the above criteria.

New features of Vicon Polygon

For information about the new features that were added to Polygon 4.3.1, see the following topics:


- [Video title display on page 7](#)
- [Time Bar displays active data folder on page 8](#)
- [Other updates on page 9](#)

Video title display



When you create a movie pane, the video name is automatically displayed as the title at the top of the pane. This is particularly useful if you are displaying multiple video files, making it easy to identify each file.

Time Bar displays active data folder



Previously, events were displayed in the Time Bar, but it was not always obvious which trial they belonged to. To make it clear which events are currently displayed, the title of the active trial (the trial currently selected in the Data Bar), is now always displayed above the Time Bar.

Other updates

- To avoid scrolling when you load a large number of trials, the panes within the **Data Bar** are now resizable.
- When you choose a custom **Graph Trace Color** from the **Attributes** panel, the **Context Appearance** is automatically disabled, so that the relevant trace color in the graph is updated.
- Memory allocation has been improved to remove the limit on the number of videos you can load.

Vicon Polygon addressed issues

This section lists issues that have been addressed in the current release of Polygon.

Issue	Description
POLY-991	Export of Protected Report (ie with Protect Exported Report selected), followed by Open Report now correctly opens the exported report.
POLY-1034	A warning message no longer appears when closing VAULT after declining to use a trial license.
POLY-1046	Searching in ProEclipse inside Polygon now correctly handles the use of a space in the search terms.
POLY-1054	Multiple pop-up error messages no longer appear on saving/closing Polygon files.

Vicon Polygon known issues

This section lists known issues in the current release of Polygon.

Issue	Solution/Workaround
Under Windows 7, if you set the text size to 125% or 150%, Polygon's version number is cropped from the About Polygon dialog box. (POLY-731).	To display the version number, change the text size. This setting can be found in Windows Control Panel > Appearance and Personalization > Display .

Contact Vicon

If you need more information than that supplied in the documentation or on the [Vicon Support web pages](#), please contact Vicon:

Denver, CO

Vicon Denver
7388 S. Revere Parkway, Suite 901
Centennial
CO 80112
USA
T:303.799.8686
F:303.799.8690
E: support@vicon.com

Oxford, UK

Vicon Oxford
14 Minns Business Park
West Way
Oxford
OX2 0JB
UK
T:+44.1865.261800
F:+44.1865.240527
E: support@vicon.com

Los Angeles, CA

Vicon LA
3750 S. Robertson Boulevard, Suite 100
Culver City
Los Angeles
CA 90066
USA
T:303.799.8686
F:310.388.3200
E: support@vicon.com

Singapore

Vicon Singapore
T:+65 6400 3500
E: support@vicon.com

© Copyright 2016 Vicon Motion Systems. All rights reserved.

