

# Modeling with MATLAB



## Important

Vicon Nexus 2 is supported with MATLAB 2013b and later. Nexus may function with other versions of MATLAB, but other versions have not been extensively tested by Vicon. To use MATLAB with Vicon Nexus 2, ensure that, in addition to installing MATLAB, you install [.NET Framework version 4.5](#).

The MATLAB interface provides immediate feedback of scripting changes in the **3D Perspective** view.

For more information, see:

- [Example of modeling with MATLAB](#)
- [Advanced MATLAB modeling](#)
- [MATLAB commands for use with Nexus](#)
- [MATLAB troubleshooting](#)