

# Vicon Nexus regulatory information

The following topics contain regulatory information relating to Vicon Nexus. For full regulatory information relating to your Vicon hardware, see the documentation that was supplied with your Vicon system.

- [ISO certification](#)
- [Medical devices directive](#)



## Important

- Vicon IMUs are supported for use with Nexus for research purposes only. For full sensor safety and regulatory details, see the [Vicon Blue Trident Model V2 Safety and Regulatory Information](#).
- The firmware for Vicon IMUs (ie, Blue Trident sensors) that are to be used with Nexus 2.10 must be updated to version 9.0.2 or later, using Capture.U Desktop 1.1, which is available from the [Capture.U page](#) on the Vicon website.

## Global Unique Device Identification Database <https://accessgudid.nlm.nih.gov/>

- Vicon Vantage systems: [B548VICONVANTAGESYSTEM0](#)
- Vicon Vero systems: [B548VICONVEROSYSTEM0](#)
- Vicon MX T-Series systems: [B548VICONTSERIESSYSTEM0](#)
- Vicon Bonita systems: [B548BONITASYSTEM0](#)