

Modeling with MATLAB



Important

Vicon Nexus 2 is compatible with, and has been tested with MATLAB R2013b. Nexus may function with other versions of MATLAB, but other versions have not been extensively tested by Vicon. To use MATLAB with Vicon Nexus 2, ensure that, in addition to installing MATLAB, you install [.NET Framework version 4.5](#)

The MATLAB interface provides immediate feedback of scripting changes in the **3D Perspective** view.

For more information, see:

- [Example of modeling with MATLAB](#)
- [Advanced MATLAB modeling](#)
- [MATLAB commands for use with Nexus](#)
- [MATLAB troubleshooting](#)